Determine™

HIV-1/2

REF 7D2346, 7D2347

Key to symbols used / Erläuterung der verwendeten Symbole /
Clave de los símbolos utilizados / Légende des symboles utilisés /
Legenda dei simboli usati / Legenda dos símbolos utilizados

Store at 2-30°C /
Lagerung bei 2 bis 30°C /
Almacenar a 2-30°C /
Conserver entre 2 et 30°C /
Conservare a 2-30°C /
Conservar a 2-30°C

EDTA CAPILLARY TUBES

EDTA Capillary Tubes / EDTA-Kapillarröhrchen /
Tubos capilares con EDTA / Tubes capillaires
contenant de l'EDTA / Microprovette con EDTA /
Tubos capilares EDTA

CHASE BUFFER

Chase Buffer / Antriebspuffer / Tampón de
arrastre / Tampon de fixation / Tampone chase /
Tampão de fixação

Alere Ltd.
Pepper Road, Hazel Grove,
Stockport, SK7 5BW, UK

0088

Alere Medical Co., Ltd.
357 Matsuhidai, Matsudo-shi,
Chiba, 270-2214 Japan
www.alere.com www.determinetest.com Tel +81 47 311 5750

EN

This package insert must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of assay results cannot be guaranteed if there are deviations from the instructions in this package insert.

INTENDED USE

For use by healthcare professionals only.
Alere Determine™ HIV-1/2 is an *In Vitro*, visually read, qualitative immunoassay for the detection of antibodies to HIV-1 and HIV-2 in human serum, plasma or whole blood. The test is intended as rapid test to detect antibodies to HIV-1/HIV-2 from infected individuals. Note: National regulations may preclude the use of rapid tests for blood screening and/or routine diagnostic analysis.

SUMMARY AND EXPLANATION OF THE TEST

AIDS (Acquired Immunodeficiency Syndrome) is characterized by changes in the population of T-cell lymphocytes. In an infected individual, the virus causes depletion of helper T-cells, which leaves the person susceptible to opportunistic infections and some malignancies. The virus that causes AIDS exists as two related types known as HIV-1 and HIV-2. The presence of the AIDS virus elicits the production of specific antibodies to either HIV-1 or HIV-2.^{1,2}

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

Alere Determine™ HIV-1/2 is an immunochromatographic test for the qualitative detection of antibodies to HIV-1 and HIV-2. Sample is added to the sample pad. As the sample migrates through the conjugate pad, it reconstitutes and mixes with the selenium colloid-antigen conjugate. This mixture continues to migrate through the solid phase to the immobilized recombinant antigens and synthetic peptides at the patient window site.
If antibodies to HIV-1 and/or HIV-2 are present in the sample, the antibodies bind to the antigen-selenium colloid and to the antigen at the patient window, forming a red line at the patient window site.
If antibodies to HIV-1 and/or HIV-2 are absent, the antigen-selenium colloid flows past the patient window, and no red line is formed at the patient window site.
To ensure assay validity, a procedural control bar is incorporated in the assay device.

CONTENTS

Alere Determine™ HIV-1/2, 20 Tests (7D2346) or 100 Tests (7D2347)

- Alere Determine™ HIV-1/2 Test Card**, 2 cards or 10 cards (10 tests/card), HIV-1/2 recombinant antigen and synthetic peptide coated.
- For testing whole blood specimens the following reagent is required:
Chase Buffer (7D2243)
- [CHASE BUFFER]** 1 Bottle (2.5 mL) Chase Buffer prepared in phosphate buffer. Preservatives: Antimicrobial Agents.

ACCESSORIES (required but not provided)

Whole Blood (fingerstick assay)

- [EDTA CAPILLARY TUBES]** (7D2227)

WARNINGS AND PRECAUTIONS

For professional use only.

CAUTION

This product contains human sourced infectious and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, it is recommended that all human sourced materials be considered potentially infectious and handled with appropriate biosafety practices.

STORAGE

Alere Determine™ HIV-1/2 Test and the Chase Buffer must be stored at 2-30°C until expiration date.
Kit components are stable until expiration date when handled and stored as directed. Do not use kit components beyond expiration date.
Immediately reseal all unused tests in the foil pouch containing the desiccant by pressing seal from end to end to close.
Do not use devices that have become wet or the packaging has become damaged.

SPECIMEN COLLECTION




Serum, Plasma and Whole Blood Collection by Venipuncture

Human serum, plasma and whole blood collected by venipuncture should be collected aseptically in such a way as to avoid hemolysis. Specimens showing particulate matter or turbidity should be centrifuged before testing in order to avoid providing erroneous results.
NOTE: For whole blood and plasma specimens, EDTA collection tubes must be used.

Whole Blood Collection by Fingertick.⁴

Before collecting a fingerstick specimen, place an EDTA capillary tube on a clean dry surface.

- Choose the fingertip of the middle, ring, or index finger (whichever is the least callused) for adults and children older than one year. Warm the hand as needed with a warm, moist towel or warm water to increase blood flow.
- Clean fingertip with alcohol; allow to air dry. Position the hand palm-side up.
- Use a new lancet for each person. Place the lancet off-center on the fingertip. Firmly press the lancet against the finger and puncture the skin. Dispose of the lancet in an appropriate biohazard sharps container.
- Wipe away the first drop of blood with a sterile gauze pad.
- Hold the finger lower than the elbow and apply gentle, intermittent pressure to the base of the punctured finger several times. Touch the tip of the EDTA Capillary Tube to the drop of blood*. Avoid air bubbles.

* If EDTA Capillary Tubes will be used, fill the tube with blood between the 2 marked lines.

SPECIMEN STORAGE

- Serum and plasma specimens should be stored at 2-8°C if the test is to be run within 7 days of collection. If testing is delayed more than 7 days, the specimen should be frozen (-20°C or colder).
- Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 7 days of collection. Do not freeze whole blood specimens.
- Whole blood collected by fingerstick should be tested immediately.

TEST PROCEDURE

The desired number of test units from the 10-test card can be removed by bending and tearing at the perforation.

NOTE

- Removal of the test units should start from the right side of the test card to preserve the lot number which appears on the left side of the test card.**
- Assay should be initiated within 2 hours after removing the protective foil cover from each test.**
- 1. Remove the protective foil cover from each test.
- 2. For serum or plasma samples:
 - Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
 - Wait a minimum of 15 minutes (up to 60 minutes) and read result.
- 3. For whole blood (venipuncture) samples:
 - Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
 - Wait one minute, then apply one drop of Chase Buffer to the sample pad.
 - Wait a minimum of 15 minutes (up to 60 minutes) and read result.
- 4. For whole blood (fingerstick) samples:
 - Apply 50 µL of sample (by EDTA capillary tube) to the sample pad (marked by the arrow symbol).
 - Wait until blood is absorbed into the sample pad, then apply one drop of Chase Buffer to the sample pad.
 - Wait a minimum of 15 minutes (up to 60 minutes) and read result.


QUALITY CONTROL

- To ensure assay validity, a procedural control is incorporated in the device and is labeled "Control". If the control bar does not turn red by assay completion, the test result is invalid and the sample should be retested.
- National regulations may require the use of both positive and negative controls with each assay.

INTERPRETATION OF RESULTS


POSITIVE (Two Bars)

Red bars appear in both the control window (labeled "Control") and the patient window (labeled "Patient") of the strip. Any visible red color in the patient window should be interpreted as positive.


NEGATIVE (One Bar)

One red bar appears in the control window of the strip (labeled "Control"), and no red bar appears in the patient window of the strip (labeled "Patient").

INVALID (No Bar)

If there is no red bar in the control window of the strip, and even if a red bar appears in the patient window of the strip, the result is invalid and should be repeated.

NOTES

- The test result is positive even if the patient bar appears lighter or darker than the control bar.
- If an invalid test result occurs repeatedly, or for technical assistance, contact your local distributor or call Technical Support.

LIMITATIONS OF THE PROCEDURE

- Alere Determine™ HIV-1/2** is designed to detect antibodies to HIV-1 and HIV-2 in human serum, plasma and whole blood. Other body fluids or pooled specimens may not give accurate results.
- The intensity of the patient bar does not necessarily correlate to the titer of antibody in the specimen.
- Because of possible non-specific reactions due to causes other than HIV infection, particularly when testing low prevalence populations, it is appropriate to further investigate specimens found to be positive. Some supplemental tests are required e.g., Western Blot (WB), Immunofluorescence Assay (IFA), Radio Immune Precipitation Assay (RIPA), Nucleic Acid Amplification Test (NAT).
- The whole blood application could not be evaluated on seroconversion panels due to the unavailability of such panels made up of whole blood specimens.
- A negative result with **Alere Determine™ HIV-1/2** does not exclude the possibility of infection with HIV.
- A false negative result can occur in the following circumstances:
 - low levels of antibody (e.g., early seroconversion specimens) are below the detection limit of the test
 - infection with a variant of the virus that is less detectable by the **Alere Determine™ HIV-1/2** assay configuration
 - HIV antibodies in the patient that do not react with specific antigens utilized in the assay configuration - in exceptional cases this may lead to observation of false negative results.
 - specimen handling conditions which result in loss of HIV antibody multivalency

For these reasons care should be taken in interpreting negative results. Other clinical data (e.g., symptoms or risk factors) should be used in conjunction with the test results.

- Positive specimens should be retested using another method and the results should be evaluated in light of the overall clinical evaluation before a diagnosis is made.
- Whole blood or plasma specimens containing anticoagulants other than EDTA may give incorrect results.

PERFORMANCE CHARACTERISTICS

SPECIFICITY

Table I Specificity of Alere Determine™ HIV-1/2		
Population	Number of Specimens Tested	Negative by Alere Determine™ HIV-1/2
Seronegative Serum	1707*	1703/1707 (99.77%)
Plasma	521	521/521 (100.00%)
Pregnant Females	208*	206/207 (99.52%)
West Africans	49	48/49 (97.96%)
Disease States Other than HIV and Potentially Interfering Substances	320*	313/316 (99.05%)
Total	2805**	2791/2800 (99.68%)

* One specimen from pregnant female, and four specimens from disease states were positive by **Alere Determine™ HIV-1/2**. Four specimens confirmed positive by Western Blot and one specimen was indeterminate by Western Blot, and positive by a commercially available test of Microparticle Enzyme Immunoassay (MEIA) method.

** 961 specimens were from Europe, 706 specimens were from North America, 1089 specimens were from Asia, and 49 specimens were from Africa.
A total of 2805 serum and plasma specimens from Europe, Asia, West Africa, and North America were tested by **Alere Determine™ HIV-1/2** (Table I).

Table II A Comparison of Alere Determine™ HIV-1/2 Specificity by Geographic Area		
Area	Number of Specimens Tested	Negative by Alere Determine™ HIV-1/2
Europe	801	796/799 (99.62%)
North America	706	697/702 (98.29%)
Asia	1089	1089/1089 (100.00%)
Africa	2118	2079/2118 (98.16%)

A total of 4714 seronegative serum and plasma specimens from Europe, North America, Asia, and Africa were tested by **Alere Determine™ HIV-1/2** (Table II). The specimens from North America, Asia, and 801 of 961 specimens from Europe, 49 of 2118 specimens from Africa (referred to as "West Africans" in Table I) were included in Table I. Discordant specimens were confirmed negative by either Western Blot or HIV-1 PCR assays.

Table III A Comparison of Alere Determine™ HIV-1/2 Specificity in Seronegative Whole Blood and Paired Serum and Plasma Specimens		
Specimen Type	Number of Specimens Tested	Negative by Alere Determine™ HIV-1/2
Serum	416*	416/416 (100.00%)
Plasma	416*	416/416 (100.00%)
Whole Blood (venipuncture)	416*	416/416 (100.00%)
Whole Blood (fingerstick)	39	39/39 (100.00%)

* 48 specimens were from Europe, 368 specimens were from Thailand.

A total of 416 seronegative whole blood specimens from Europe and Thailand were tested with paired serum and plasma by **Alere Determine™ HIV-1/2**. Thirty-nine of the whole blood specimens from Thailand were collected by both venipuncture and fingerstick (Table III).

SENSITIVITY

Table IV Sensitivity of Alere Determine™ HIV-1/2		
Population	Number of Specimens Tested	Positive by Alere Determine™ HIV-1/2
HIV-1 Positive	599*	599/599 (100.00%)
HIV-2 Positive	147*	147/147 (100.00%)
HIV-1 Subtypes A-G	222	222/222 (100.00%)
HIV-1 Group O	12	12/12 (100.00%)
Total	980	980/980 (100.00%)

* 111 specimens were from Europe, 228 specimens were from North America, 296 specimens were from Asia, and 111 specimens were from Africa.

A total of 980 HIV-1 and HIV-2 antibody positive serum and plasma specimens from Europe, Asia, Africa, North and South America were tested by **Alere Determine™ HIV-1/2** (Table IV).

Table V A Comparison of Alere Determine™ HIV-1/2 Sensitivity by Geographic Area		
Area	Number of Specimens Tested	Positive by Alere Determine™ HIV-1/2
Europe	111	111/111 (100.00%)
North America	228	228/228 (100.00%)
Asia	296	296/296 (100.00%)
Africa	1129	1128/1129 (99.91%)

* One negative specimen by **Alere Determine™ HIV-1/2** confirmed positive by HIV-1 PCR and a commercially available test of MEIA method.

A total of 1764 seropositive serum and plasma specimens from Europe, North America, Asia, and Africa were tested by **Alere Determine™ HIV-1/2** (Table V). The specimens from Europe, North America, Asia, and 111 of 1129 specimens from Africa (referred to as "HIV-2 Positive" in Table IV) were included in Table IV. Discordant specimens were confirmed HIV-1 positive by either Western Blot or HIV-1 PCR assays.

Table VI A Comparison of Alere Determine™ HIV-1/2 Sensitivity in Seropositive Whole Blood and Paired Serum and Plasma Specimens		
Specimen Type	Number of Specimens Tested	Positive by Alere Determine™ HIV-1/2
Serum	127*	127/127 (100.00%)
Plasma	127*	127/127 (100.00%)
Whole Blood (venipuncture)	127*	127/127 (100.00%)
Whole Blood (fingerstick)	32	32/32 (100.00%)

* 25 specimens were from Europe, 102 specimens were from Thailand.

A total of 127 seropositive whole blood specimens from Europe and Thailand were tested with paired serum and plasma by **Alere Determine™ HIV-1/2**. Thirty-two of the whole blood specimens from Thailand were collected by both venipuncture and fingerstick (Table VI).

Table VII A Comparison of Alere Determine™ HIV-1/2 Sensitivity in Seroconversion panel with a Commercially Available Test			
Panels	Day of Bleed	Commercially Available Test* (S/CO)	Alere Determine™ HIV-1/2
Panel A			
A-4	18	0.55	-
A-5	21	1.45	-
A-6	25	12.84	+
A-7	28	16.85	+
Panel B			
B-4	9	0.38	-
B-5	15	0.38	-
B-6	28	6.59	+
B-7	33	12.18	+
Panel C			
C-3	13	0.39	-
C-4	27	1.47	-
C-5	34	8.76	+
C-6	50	5.32	+
Panel D			
D-1	0	0.38	-
D-2	3	0.39	-
D-3	9	5.76	+
Panel E			
E-1	0	0.38	-
E-2	7	0.54	-
E-3	11	4.59	+
E-4	15	6.09	+

* The reference method of a commercially available test is MEIA.

5 seroconversion panels were tested by **Alere Determine™ HIV-1/2** and a commercially available test of MEIA method (Table VII).

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