Key to symbols used / Schlüssel für die verwendeten Symbole /  
Clave de los símbolos utilizados / Légende des symboles utilisés /  
Legenda dei simboli utilizzati / Designação dos símbolos utilizados

REF

Catalogue Number / Katalognummer /  
Número de catálogo / Référence  
catalogue / Numero di catalogo /  
N° de Catálogo

IVD

*In Vitro* Diagnostic Medical Device /  
*In Vitro* - Diagnostika / Dispositivo  
de diagnóstico médico *in vitro* /  
Dispositif médical de diagnostic  
*In Vitro* / Dispositivo medico  
diagnostico *In Vitro* / Dispositivo  
Médico para Diagnóstico *In Vitro*

2°C

30°C

Store at 2-30°C / Bei 2-30°C lagern /  
Guardar a temperaturas entre 2 y  
30°C / Conserver entre 2-30°C /  
Conservare a 2-30°C / Armazenar  
a 2-30°C

Σ

20

Contains Sufficient for 20 tests /  
Enthält eine ausreichende Menge  
für 20 Tests / Contiene material  
suficiente para realizar 20 pruebas /  
Permet de réaliser 20 tests /  
Contiene reagenti sufficienti per 20  
tests / Contém o suficiente para 20  
testes

Σ

100

Contains Sufficient for 100 tests /  
Enthält eine ausreichende Menge  
für 100 Tests / Contiene material  
suficiente para realizar 100 pruebas /  
Permet de réaliser 100 tests /  
Contiene reagenti sufficienti per 100  
tests / Contém o suficiente para 100  
testes

Alere Medical Co., Ltd.  
357 Matsuhidai, Matsudo-shi,  
Chiba, 270-2214 Japan  
[www.alere.com](http://www.alere.com) [www.determinetest.com](http://www.determinetest.com) Tel +81 47 311 5750

February 2012  
240795/R6

**LIMITATIONS OF THE PROCEDURE**

- The **Alere** Determine™ HIV - 1/2 test is designed to detect antibodies to HIV - 1 and HIV - 2 in human serum, plasma, and whole blood. Other body fluids or pooled specimens may not give accurate results.
- The intensity of the patient bar does not necessarily correlate to the titer of antibody in the specimen.
- A negative result with **Alere** Determine™ HIV-1/2 does not exclude the possibility of infection with HIV. A false negative result can occur in the following circumstances:
  - low levels of antibody (e.g., early seroconversion specimens) are below the detection limit of the test
  - infection with a variant of the virus that is less detectable by the **Alere** Determine™ HIV-1/2 assay configuration
  - HIV antibodies in the patient that do not react with specific antigens utilized in the assay configuration
  - specimen handling conditions which result in loss of HIV antibody multivalency

For these reasons care should be taken in interpreting negative results. Other clinical data (e.g., symptoms or risk factors) should be used in conjunction with the test results.

- Positive specimens should be retested using another method and the results should be evaluated in light of the overall clinical evaluation before a diagnosis is made.
- Whole blood or plasma specimens containing anticoagulants other than EDTA may give incorrect results.

PERFORMANCE CHARACTERISTICS

SPECIFICITY

A total of 1594 serum and plasma specimens from Asia, West Africa, and North America were tested by **Alere** Determine™ HIV-1/2 and a commercially available test (Table I).

Table I Specificity of <b>Alere</b> Determine™ HIV-1/2			
Population	Number of Specimens Tested	Negative by <b>Alere</b> Determine™ HIV-1/2	Negative by a Commercially Available Test***
Seronegative Serum	908	907/908 (99.89%)	908/908 (100.00%)
Plasma	403	403/403 (100.00%)	403/403 (100.00%)
Pregnant Females	58*	57/57 (100.00%)	57/57 (100.00%)
West Africans	49	48/49 (97.96%)	48/49 (97.96%)
Disease States Other than HIV and Potentially Interfering Substances	176*	173/175 (98.86%)	174/175 (99.45%)
Total	1594**	1588/1592 (99.75%)	1590/1592 (99.87%)

\* One specimen from a pregnant female and an HCV positive patient were positive by both **Alere** Determine™ HIV-1/2 and the commercially available test. Both specimens confirmed positive by HIV-1 Western Blot.

\*\* 456 specimens were from North America, 1089 specimens were from Asia, and 49 specimens were from Africa.

\*\*\* The reference method of a commercially available test is particle agglutination.

WARNINGS AND PRECAUTIONS

For *In Vitro* Diagnostic Use.

**CAUTION:** Appropriate biosafety practices<sup>4,5</sup> should be used when handling specimens and reagents. These precautions include, but are not limited to the following:

- Wear gloves.
- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect all spills of specimens or reagents using a suitable disinfectant, such as 0.5% sodium hypochlorite.<sup>6,7</sup>
- Decontaminate and dispose of all specimens, reagents, and other potentially contaminated materials in accordance with local regulations.<sup>8,9</sup>

STORAGE

The **Alere** Determine™ HIV-1/2 Test Cards and Chase Buffer must be stored at 2-30°C until expiration date.

Kit components are stable until expiration date when handled and stored as directed. Do not use kit components beyond expiration date.

SPECIMEN COLLECTION

Serum, Plasma, and Whole Blood Collection by Venipuncture

Human serum, plasma, and whole blood collected by venipuncture should be collected aseptically in such a way as to avoid hemolysis.

**NOTE: For whole blood and plasma specimens, EDTA collection tubes must be used.**

**Whole Blood Collection by Fingertick<sup>10</sup>**

Before collecting a fingertick specimen, place an EDTA capillary tube on a clean dry surface.


1. Choose the fingertip of the middle, ring, or index finger (whichever is the least callused) for adults and children older than one year. Warm the hand as needed with a warm, moist towel or warm water to increase blood flow.

2. Clean fingertip with alcohol; allow to air dry. Position the hand palm-side up.

3. Use a new lancet for each person. Place the lancet off-center on the fingertip. Firmly press the lancet against the finger and puncture the skin. Dispose of the lancet in an appropriate biohazard sharps container.

4. Wipe away the first drop of blood with a sterile gauze pad.

5. Hold the finger lower than the elbow and apply gentle, intermittent pressure to the base of the punctured finger several times. Touch the tip of the EDTA Capillary Tube to the drop of blood\*. Avoid air bubbles.

\*If EDTA Capillary Tubes (No. 7D2222) will be used, fill the tube with blood between the 2 marked lines.

SPECIMEN STORAGE

- Serum and plasma specimens should be stored at 2-8°C if the test is to be run within 7 days of collection. If testing is delayed more than 7 days, the specimen should be frozen (-20°C or colder).
- Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 7 days of collection. Do not freeze whole blood specimens.
- Whole blood collected by fingertick should be tested immediately.

TEST PROCEDURE

The desired number of test units from the 10-test card can be removed by bending and tearing at the perforation.

**NOTE: Removal of the test units should start from the right side of the test card to preserve the lot number which appears on the left side of the test card.**

1. Remove the protective foil cover from each test.
2. For serum or plasma samples:
  - a. Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
  - b. Wait a minimum of 15 minutes (up to 60 minutes) and read result.
3. For whole blood (venipuncture) samples:
  - a. Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
  - b. Wait one minute, then apply one drop of Chase Buffer to the sample pad.
  - c. Wait a minimum of 15 minutes (up to 60 minutes) and read result.
4. For whole blood (fingertick) samples:
  - a. Apply 50 µL of sample (by EDTA capillary tube) to the sample pad (marked by the arrow symbol).
  - b. Wait until blood is absorbed into the sample pad, then apply one drop of Chase Buffer to the sample pad.
  - c. Wait a minimum of 15 minutes (up to 60 minutes) and read result.

QUALITY CONTROL

To insure assay validity, a procedural control is incorporated in the device and is labeled "Control". If the control bar does not turn red by assay completion, the test result is invalid and the sample should be retested.

INTERPRETATION OF RESULTS

POSITIVE (Two Bars)

Red bars appear in both the control window (labeled "Control") and the patient window (labeled "Patient") of the strip. Any visible red color in the patient window should be interpreted as positive.


NEGATIVE (One Bar)

One red bar appears in the control window of the strip (labeled "Control"), and no red bar appears in the patient window of the strip (labeled "Patient").

INVALID (No Bar)

If there is no red bar in the control window of the strip, and even if a red bar appears in the patient window of the strip, the result is invalid and should be repeated.

NOTES:

- The test result is positive even if the patient bar appears lighter or darker than the control bar.
- If an invalid test result occurs repeatedly, or for technical assistance, contact your local distributor or call Technical Support.

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

1. Plot P, Plummer FA, Mhalu FS, Lamborary JL, Chin J, Mann JM. AIDS: An International Perspective. Science. 1988;239:573-579.
2. Weniger BG, Takebe Y, Ou CY, et al. The Molecular Epidemiology of HIV in Asia. AIDS. 1994;8(S2):S13-S28.
3. Gürtler LG, Hauser PH, Eberle J, et al. A New Subtype of Human Immunodeficiency Virus Type 1 (MVP-5180) from Cameroon. Journal of Virology. 1994;68(3):1581-1585.
4. World Health Organization. Laboratory Biosafety Manual. Geneva: World Health Organization,1993.
5. National Committee for Clinical Laboratory Standards. Protection of Laboratory Workers from Infectious Disease Transmitted by Blood, Body Fluids, and Tissue. Tentative Guideline. NCCLS Document M29-T2. Villanova, PA: NCCLS, 1991:1-43.
6. CDC. Recommendations for Prevention of HIV Transmission in Health-Care Settings, MMWR 1987;36(2S):3S-18S.
7. Sehustler LM, Hollinger FB, Dreesman GR, et al. Immunological and Biophysical Alteration of Hepatitis B Virus Antigens by Sodium Hypochlorite Disinfection. Appl Envir Microbiol 1981;42(5):762-7.
8. National Committee for Clinical Laboratory Standards. Clinical Laboratory Waste Management: Approved Guideline. NCCLS Document GP5-A. Villanova, PA: NCCLS, 1993;13(22):1-18, 29-42.
9. US Environmental Protection Agency. EPA Guide for Infectious Waste Management. Publication No. EPA/530-SW-86-014. Washington, DC: US Environmental Protection Agency, 1986:1-1-5-5, R1-R3, A1-A24.
10. NCCLS Guidelines, No. H4A3, Volume 11, No. 11, Approved Standard, July 1991.

Advice Line

For further information, please contact your distributor, or call to one of the following Alere Product Support Care Centers:

Region	Phone	E-Mail Address
Europe & Middle East	+ (44) 161 483 9032	EMEproductsupport@alere.com
Asia Pacific	+ (61) 7 3363 7711	APproductsupport@alere.com
Africa, Russia & CIS	+ (972) 8 9429 683	ARCISproductsupport@alere.com
Latin America	+ (57) 01800 094 9393	LAprодуктsupport@alere.com

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